

## Sophisticated Positioning Aids for conventional X-Ray



⊕ Perfectly hygienic

⊕ High durability

⊕ Versatile Precision

### What our customers say

“Even in the most difficult situations: the positioning aids **remain intact** and enable **versatile** and **precise** imaging.

In **unsanitary conditions** I no longer have to **wrap** my sponges and can even position high-risk patients without any **hygiene concerns**. Patient positioning is **fun** with these tools.”

- Michael Bunte,  
X-Ray Advisor, Germany

### Our Solution

The positioning aids are made of a hard-wearing PU foil and are sealed using a high-frequency welding process. As a result, they withstand the high stresses of the radiology department and can be cleaned even with strong disinfectants.





### ProFoam

Comprehensive range with various standard shapes and sizes:

- Smooth, easily disinfectable surfaces
- Welded seams
- Optional: extra-stable foam material for precise imaging with heavy patients

### ProFoam Table Mat

ProFoam Table Mats are available in different dimensions and designs:

- Slideable or anti-slip bottom side
- With or without handles



### Leg Rest

Professional and comfortable elevated positioning with leg support stand made of acrylic glas, for e.g. cross table lateral hip projection:

- Height-adjustable
- Molded support shell

### Further Applications

