

Product Data Sheet



Sandbag Positioning Aid

Ref.:	Product name:	Dimensions:	Weight:	Packaging dimensions:
1261-A	Sandbag 40x20, 3kg	400 x 200 x 50 mm	3000 g	280 x 470 mm
1262-A	Sandbag 30x20, 2kg	300 x 200 x 50 mm	2000 g	280 x 470 mm
1263-A	Sandbag 24x18, 1.25kg	240 x 180 x 40 mm	1250 g	280 x 470 mm

Intended use	Short-term patient positioning and stabilization device during invasive and non-invasive procedures to be used on intact patient's skin. Non-Sterile and not to be used in the sterile field.		
Directives & classification	CE according to MDR – Regulation 2017/745 Class 1 product		
Applied quality assurance	SN EN ISO 13485:2016		
Instructions for use	Enclosed to the product		
Restrictions of use	Do not use the product for purposes other than those described.		
General safety notes	Non sterile. Do not sterilize. Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).		
MR safety note	MR Safe		
Materials	Material TPU Quartz sand LDPE	Component foils sand packaging	
Marking	CE		
Storage	Store in dry place / room temperature		

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.